


# Measurement of the Potential Exposure to Isocyanate Monomers and Oligomers During Spray Painting

**Gary Y Oishi, Jamie L. Brown, Emily R. Barrey, Olga I. Shimelis, Michael R. Halpenny, Kristen L. Schultz**  
Supelco/Sigma-Aldrich, 595 North Harrison Road, Bellefonte PA 16823




[sigma-aldrich.com/analytical](http://sigma-aldrich.com/analytical)

# What is an isocyanate? (chemistry)



*Isocyanate + polyol (compound with many alcohol groups) = polyurethane*



## Isocyanates In the News:

January 8, 2015 – US EPA issues rule proposal limiting the US of TDI chemicals in consumer products.

2014 – Worker exposure to TDI at Automotive manufacturing plant in Selma, LA gains national attention

June 25, 2013 – OSHA NEP on isocyanate exposure in the workplace

April 2013 – 1,2-MP and other piperazine banned in UK as Class C drugs of abuse.

April 13, 2011 – EPA Action Plan on MDI & TDI



# Isocyanates Exposure

Isocyanates exposure is universal

## Health Risks

- Respiratory: Asthma - **leading cause of occupational asthma**; difficulty breathing [Johnson et al. 2004; Wisnewski et al. 2006]
- Musculoskeletal: Chest pains
- Irritation of mucous membranes; dermatitis
- Cancer: Aromatic isocyanates (Europe)
- Death: MIC, Union Carbide (Bhopal, India)

## NIOSH Pocket Guide to Chemical Hazards

Hexamethylene diisocyanate (HDI) (CAS No. 822-06-0)

Isophorone diisocyanate (IPDI) (CAS No. 4098-71-9)

Methyl isocyanate (MIC) (CAS No. 624-83-9)

Methylenebis(phenyl isocyanate)(MDI) (CAS No. 101-68-8)

Naphthalene diisocyanate (NDI) (CAS No. 3173-72-6)

Toluene-2,4-diisocyanate (TDI) (CAS No. 584-84-9)

| Isocyanate Species                     | NIOSH REL (ug/m3)    | ACGIH TLV (ug/m3) |
|--|----------------------|-------------------|
| TDI                                    | CA (lowest feasible) | 36                |
| MDI                                    | 50                   | 51                |
| HDI                                    | 35                   | 34                |
| HMDI                                   | None                 | 54                |
| IPDI                                   | 45                   | 45                |
| NDI                                    | 40                   | None              |
| TRIG (total reactive isocyanate group) | None                 | None              |

# Facts

Isocyanates are a main component in the production of polyurethane (PUR) materials

Exposure to isocyanates can put workers at risk for respiratory disorders like “occupational asthma”

Not only monomers pose a health risk

# How Isocyanate Exposure Occurs

Application – before curing is complete



Thermal degradation – Grinding, Cutting, Welding





## Current Test Methods

### NIOSH Method 2535 – TDI

- Glass tube with (N-[(4-nitrophenyl)methyl]-propylamine on glass wool
- HPLC UV

### ASTM D5932 – TDI

- Segregating device containing with a glass fiber filter impregnated with 9-(N-methylaminomethyl) anthracene
- HPLC UV/Fluorescence

### ASTM D6561 – oligomeric HDI

- Segregating device containing a polytetrafluoroethylene (PTFE) filter
- HPLC UV

### OSHA 42/47 - TDI, HDI, MDI

- Filter cassette with glass fiber filter impregnated with 1-(2-pyridyl)piperazine
- HPLC UV

### ISO 17734-1:2013 – monomeric and oligomeric

- Filter cassette + denuder with glass fiber filter impregnated with dibutylamine
- LC-MS/MS





## Experiment-Overview

Spray painting using an automotive clear coat

Painting was done for a continuous length of time

Collect personal and area samples

Quantitative analysis was performed accordingly to the specific sampler



## Spray Booth Samples

Two part clear coat was used to spray paint a car hood





## Spray Booth Samples

Two different types of “dry” samplers



### ORBO-80

1-(2-pyridyl)piperazine (1-2PP) coated glass fiber filter  
1L/min per OSHA 42/47



### Asset EZ4-NCO

Dibutylamine (DBA) coated glass fiber filters  
100ml/min



# Spray Booth Samples

Two types of samples (Personal and Area)






## Spray Booth Samples

The sampling tree was placed in the rear of the spray booth



### Rear of Spray Booth



Not to scale




## Spray Booth Samples

Sampling rates:


- ORBO-80 – 1 L/min per OSHA 42/47
- Asset EZ4-NCO – 100 ml/min

- **Sampling Time: 15 minutes**






# Chromatogram – ORBO-80



LC-UV chromatogram of the resulting extract from the ORBO-80 sampler used for spray booth sampling.



# Chromatogram - ASSET



LC-MS/MS chromatogram for an ASSET area sampler used for spray booth sampling





# Calibration Standards

## ASSET - CRM

### Monomer: Dibutylamine derivative

- EIC (Ethyl isocyanate)
- HDI (Hexamethylene diisocyanate)
- HMDI (Dicyclohexylmethane 4,4'-diisocyanate)
- ICA (Isocyanic acid)
- IPDI (Isophorone diisocyanate)
- MIC (Methyl isocyanate)
- 4,4' -MDI (Methylenediphenyl diisocyanate)
- PhI (Phenyl isocyanate)
- PIC (Propyl isocyanate)
- 2,4 and 2,6 -TDI (Toluene diisocyanate)
- m-TMXDI (m-tetramethylxylene diisocyanate)

### Oligomers: Dibutylamine derivative

- HDI-Biuret
- HDI-Isocyanurate
- HDI-Uretdione
- IPDI-Isocyanurate
- MDI-Trimer
- MDI-Tetramers

## ORBO-80


### Monomers: 1-(2-pyridyl)piperazine derivative

- HDI (Hexamethylene diisocyanate)
- 4,4' -MDI (Methylenediphenyl diisocyanate)
- 2,4 and 2,6 -TDI (Toluene diisocyanate)





# Results from Personal Samples



HDI ND ( $< 4.5 \mu\text{g}/\text{m}^3$ )

HDI-Isocyanurate  $1081 \mu\text{g}/\text{m}^3$   
IPDI-Isocyanurate  $248 \mu\text{g}/\text{m}^3$


HDI-Biuret  $4 \mu\text{g}/\text{m}^3$   
HDI-Isocyanurate  $1797 \mu\text{g}/\text{m}^3$   
HDI-Uretidone  $9 \mu\text{g}/\text{m}^3$   
IPDI-Isocyanurate  $567 \mu\text{g}/\text{m}^3$

ASSET

ORBO-80



## Results from Area Samples



HDI

Conc:  $\mu\text{g}/\text{m}^3$

6.3

29.7

25.3

6.7

31.8

26.5

5.7

Red squares represent ASSET samplers, grey circles are ORBO-80 samplers.

## Results ( $\mu\text{g}/\text{m}^3$ ) from Area Samples(ASSET)


| Analyte           | M.W.  | Sample 1 | Sample 2 | Sample 3 | Sample 4 |
|-------------------|-------|----------|----------|----------|----------|
| HDI               | 168.2 | 25.3     | 29.7     | 26.5     | 31.8     |
| HDI-Biuret        | 478.6 | 21.5     | 26.9     | 22.8     | 27.1     |
| HDI-Isocyanurate  | 504.6 | 7022.4   | 8602.1   | 7796.8   | 9866.4   |
| HDI-Uretidone     | 336.4 | 50.2     | 61.3     | 51.5     | 66.4     |
| IPDI              | 222.3 | 8.6      | 8.9      | 8.5      | 10.0     |
| IPDI(2)           | 222.3 | 2.3      | 2.6      | 2.4      | 3.0      |
| IPDI-Isocyanurate | 666.9 | 5160.8   | 5801.4   | 5177.9   | 6887.9   |

## Results ( $\mu\text{g}/\text{m}^3$ ) from Area Samples(ORBO-80)


| Analyte | Sample 1 | Sample 2 | Sample 3 |
|---------|----------|----------|----------|
| 2,6-TDI | ND<4.5   | ND<4.6   | ND<4.5   |
| 1,6-HDI | 6.3      | 6.7      | 5.7      |
| 2,4-TDI | ND<4.5   | ND<4.6   | ND<4.5   |
| 4,4-MDI | ND<5.6   | ND<5.7   | ND<5.6   |

# The ASSET Sampler Works Differently

Typically the filter is placed before the vapor collection device when sampling for both vapors and particulates.




The ASSET Sampler works Differently.....the filter is placed after the vapor collection device (denuder). The vapors collect in the denuder, while the particles pass through to the filter.






Environmental

# Reagent Depletion – During Sampling



Filter



37mm Filter Cassette



Denuder

Filter

## Limit of Quantification: LC-MS/MS vs. LC-MS

*Calibration concentration range: 5-280 ng/mL*

| Compound | MS (ng/mL) | MS/MS (ng/mL) |
|----------|------------|---------------|
| ICA      | 5.00       | 1.000         |
| MIC      | 25.00      | 3.000         |
| EIC      | 15.00      | 4.000         |
| PIC      | 10.00      | 10.000        |
| PhI      | 3.00       | 0.400         |
| HDI      | 2.00       | 0.500         |
| 2,6-TDI  | 3.00       | 0.030         |
| 2,4-TDI  | 3.00       | 0.030         |
| IPDI-1   | 4.00       | 0.400         |
| IPDI-2   | 4.00       | 0.400         |
| MDI      | 3.00       | 2.000         |





## Key Differences

- Calibration standards
- Analytical technique
- Sampler design

**For More Information on the Dry Sampler, visit our website or request the flyer.**

**4-Page Flyer – T412089 (OTB)**

**On the web at: [www.sigmaaldrich.com/as](http://www.sigmaaldrich.com/as)**



For more information, visit **[sigmaaldrich.com](http://sigmaaldrich.com)**

**SUPELCO**  
Sigma-Aldrich

**ASSET™ EZ4-NCO Dry Sampler for Isocyanates**  
First Class Sampling and Analysis of Isocyanates

**SIGMA-ALDRICH**

**ASSET EZ4-NCO Dry Sampler for Isocyanates**

**In the News**

- The ASSET EZ4-NCO automated sampler is a first-class solution for the collection and analysis of isocyanates in the workplace.
- The ASSET EZ4-NCO automated sampler is a first-class solution for the collection and analysis of isocyanates in the workplace.
- The ASSET EZ4-NCO automated sampler is a first-class solution for the collection and analysis of isocyanates in the workplace.

**About the ASSET EZ4-NCO Dry Sampler**

The ASSET EZ4-NCO automated sampler is a first-class solution for the collection and analysis of isocyanates in the workplace. It is designed to collect and analyze isocyanates in the workplace, providing a safe and accurate method for sampling and analysis.

**Advantages of the ASSET Sampler**

- Simple and safe to use
- Easy to set up and use
- Accurate and reliable results
- Low maintenance and cost
- Wide range of applications

**Key Features**

- Automated sampling and analysis
- Real-time monitoring and control
- Easy to use and maintain
- Wide range of applications

**SIGMA-ALDRICH**

**Certified Reference Materials**

The Sigma-Aldrich Certified Reference Materials (CRM) are available in a wide range of quantities and formats. They are used for the calibration and validation of analytical instruments and methods.

| Material         | Quantity | Price    | Availability |
|------------------|----------|----------|--------------|
| Isocyanate CRM 1 | 100 mg   | \$100.00 | Available    |
| Isocyanate CRM 2 | 100 mg   | \$100.00 | Available    |
| Isocyanate CRM 3 | 100 mg   | \$100.00 | Available    |

**SUPELCO**  
Sigma-Aldrich

**Analytical Service Profiles**

| Service                | Details                                   | Benefits                       | Capacity                  | Notes                                       |
|------------------------|---|--------------------------------|---------------------------|---|
| Sample Analysis        | Analysis of isocyanates in the workplace. | Accurate and reliable results. | 100 samples per day.      | See Sigma-Aldrich website for more details. |
| Instrument Calibration | Calibration of analytical instruments.    | Accurate and reliable results. | 100 calibrations per day. | See Sigma-Aldrich website for more details. |
| Method Development     | Development of analytical methods.        | Accurate and reliable results. | 100 methods per day.      | See Sigma-Aldrich website for more details. |

**SIGMA-ALDRICH**

# Conclusions

- The possible improvement that sampler design could provide
- The LC-MS/MS analysis detected and confirmed the identity of isocyanate compounds and provided accurate quantitation with the use of derivative specific and deuterated derivative standards
- Predominant airborne species (oligomers) were HDI-isocyanurate and IPDI-isocyanurate
- Higher levels of oligomers in comparison to monomer levels



# Thank You !

# Questions ?

[gary.oishi@sial.com](mailto:gary.oishi@sial.com)